

## 1.10987 Phenol Red Broth Base

Powder

Powder			
			For general questions please contact our Customer Service:
			EMD Chemicals, Inc 480 S. Democrat Road 08027 Gibbstowr United States Phone: 800 854 3417 Fax: 800 776 0999 E-mail: emdinfo@emdchemicals.com
			26 July 201
Product number	Packaging	Size	Price / Availability
1.10987.0500	Poly Bottle	500 g	USD 107.06
			Prices are subject to change without notice
Product information			
Chemical formula	Mixture		
Suitability	For Microbiology		
Specifications			
Expiration date		Exp 60 months from r	nfg date
Test date		Per Merck CofA	
Appearance (clearness)		Clear	
Appearance (color)		Red	
Color change: yellow E. faecalis ATCC 11700		+	
Color change: yellow K. pneumoniae ATCC 13883		+	
Color change: yellow Proteus vulgaris ATCC 13315		+	
Color change: yellow S. aureus ATCC 25923		+	
Color change: yellow S. typhimurium ATCC 14028		-	
Color change: yellow Shigella flexneri ATCC 12022		+/- (orange)	
Gas- Enterococcus faecalis ATCC 11700		-	
Gas- Klebsiella pneumoniae ATCC 13883		+	
Gas- Proteus vulgaris ATCC 13315		+ (poor)	
Gas- Salmonella typhimurium ATCC 14028		-	
Gas- Shigella flexneri ATCC 12022		-	
Gas- Staphylococcus aureus ATCC 25923		-	
Growth- Enterococcus faecalis ATCC 11700		Good to very good	

Growth- Klebsiella pneumoniae ATCC 13883 Growth- Proteus vulgaris ATCC 13315 Growth- Salmonella typhimurium ATCC 14028 Growth- Shigella flexneri ATCC 12022 Growth- Staphylococcus aureus ATCC 25923 pH at 25°C Good to very good Good to very good Good to very good Good to very good 7.2-7.6

## Supplemental information

Typical Composition (g/L)	Peptone from casein 5.0; peptone from meat 5.0; sodium chloride 5.0; phenol red 0.018.	
Also to be added (g/L)	Reactant 5.0-10.0 (e.g., glucose, maltose, etc.); if required, agar-agar: 0.5-1.0.	
Application	Test culture medium used together with various reactants for the biochemical identification of microorganisms by means of fermentation test.	
Preparation	15.0g/L	
Suitability	For Microbiology	

© EMD Chemicals Inc. is the North American affiliate of Merck KGaA, Darmstadt, Germany, emdinfo(at)emdchemicals.com, 2010 All references to Merck refer to Merck KGaA, Darmstadt, Germany.